



NBIMCC Medium 462

MEDIUM WITH SOIL EXTRACT

Medium 129 made up with soil extract instead of water.

Soil extract: Sterilize 1.0 kg of air-dried garden soil in 1.0 l of tap water at 121°C for 30 min. Allow to sediment for a few hours at room temperature. Add 2.0 g CaCO₃, filter and make filtrate up to 1.0 l with tap water.

Sterilize at 121°C for 30 min.

NBIMCC Medium 129

Nutrient Agar IV

Yeast extract	2.0 g
Meat extract	1.0 g
Peptone	5.0 g
NaCl	5.0 g
Agar	15.0 g
Distilled water	1000 ml

pH 7.2-7.4