



NBIMCC Medium 519

MEDIUM FOR *BIFIDOBACTERIUM*

Casein peptone, tryptic digest	10.0 g
Meat extract	5.0 g
Yeast extract	5.0 g
Glucose	10.0 g
K ₂ HPO ₄	3.0 g
Tween 80	1.0 g
Distilled water	1000 ml

pH 6.8

After sterilization, aseptically add solutions of sodium ascorbate and L-cysteine·HCl to final concentration of 1.0% and 0.05%, respectively. Heat medium not freshly prepared in a steamer for 10 min before addition of the reducing substances.